

wellion[®]

MICRO-PUMP

INSULIN PUMP SYSTEM



SMALL - EASY - DISCRETE


wellion[®]

MICRO-PUMP

INSULIN PUMP SYSTEM

FREEDOM & WELL-BEING

You decide where you want to carry your pump!



The **tubeless** Wellion[®] MICRO-PUMP can be attached at many different parts of the body to match your daily schedule or outfit. (e.g. belly, upper arm, thigh, lower back, ...).

EVERYTHING AT HAND

Diabetes can be perfectly controlled with the wireless **Portable Diabetes Assistant (PDA)**. It helps meeting your insulin needs and collects important data relating to your diabetes at the same time.


With the **integrated blood glucose meter** you can determine your blood sugar level at any time and the integrated bolus calculator calculates a suggested bolus amount.



Your AUSTRIAN health brand

DISCRETE INJECTION OF INSULIN

Due to the **small and flat design**, the insulin pump can also be worn very discreetly.



Discreet insulin delivery with the **Quick Bolus button** directly on the pump.
Step size per button press: 0.1-2U



← **Quick Bolus Button**

wellion[®]

MICRO-PUMP

INSULIN PUMP SYSTEM

INNOVATIVE DIABETES CONTROL

Wellion[®] MICRO-PUMP - the best out of patch pumps and traditional insulin pumps.

The PDA impresses with its **intuitive menu navigation** and enables extensive **data management** through the tailor-made app.



RESOURCE-FRIENDLY AND PRECISE

Thanks to the **rechargeable battery** of the Wellion[®] MICRO-PUMP, only the infusion set has to be replaced at least every 72 hours. The pump itself can be used for up to 4 years.

The latest technology of the Wellion[®] MICRO-PUMP enables precise insulin dosing starting at 0.025 U.

w wellion®



Your AUSTRIAN health brand

BLOOD GLUCOSE METER ALWAYS AT HAND

The **Wellion® MICRO-PUMP** system consists of two parts:
tubeless pump and Portable Diabetes Assistant (PDA)



wellion[®]

MICRO-PUMP

INSULIN PUMP SYSTEM

ADVANTAGES AT A GLANCE

Bolus Increment:	0.025 / 0.05 / 0.5 / 0.1 / 1U
Max Bolus Size:	35U
Extended Bolus:	% or U
Basal:	7 Basal Programs 6 or 24 U
Maximum Basal Rate:	35 U/h
Temporary Basal Rate:	% or U
Bolus Calculator:	Active insulin
	48 ranges for Target blood glucose level
	48 ranges for the ratio Insulin:Carbohydrates
	48 ranges for correction factor
Alarm notification:	Sound, light, vibration
Blood glucose test strips:	Wellion MICRO test strips
Sample volume:	0,5 µL
Measuring time:	5 seconds



PDA



MICRO-PUMP



Reservoir



Pump base

Cannula




Your AUSTRIAN health brand

CONSUMABLES



INSULIN RESERVOIR

The reservoir packaging includes:



The maximum capacity of the reservoir is **2,0 ml** → 200 units of U-100 insulin



PUMP BASE

- disposable items
- stretchy
- breathable



CANNULA

- 90° soft cannula
- 6mm / 9mm length
- unique design



Infusion site:

- belly
- upper arm
- inner thighs
- lower back
- hip

The tubeless pump is small, very lightweight and can **deliver insulin continuously for up to 3 days**, then the cannula and, if necessary, the reservoir must be changed.

Patients can shape their lives individually:
when moving, sleeping, playing or doing sports

START EASILY



1. fill the reservoir



2. attach the
infusion set



3. apply the pump



4. control with PDA

wellion[®]

MICRO-PUMP

INSULIN PUMP SYSTEM

AUSTRIA
Your
AUSTRIAN
health brand

THE BEST OUT OF PATCH PUMPS AND TRADITIONAL INSULIN PUMPS



discrete



72 hours insulin delivery



tubeless



high wearing comfort



data management



flexible lifestyle



user-friendly PDA



safe and reliable



less injections

Service and Distribution:



MED TRUST Handelsges.m.b.H.

Gewerbepark 10

7221 Marz, AUSTRIA

www.medtrust.at & www.wellion.at

Tel.: +43 2626 64 190

Fax: +43 2626 64 190 - 77

E-Mail: office@medtrust.at